

GCM2 antibody

Product Information

Catalog No.:	FNab03388
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Probable transcriptional regulator.

Immunogen information

Immunogen:	glial cells missing homolog 2(Drosophila)
Synonyms:	Chorion-specific transcription factor GCMb (hGCMb) GCM motif protein 2 Glial cells missing homolog 2 GCM2 GCMB
Observed MW:	50-56 kDa
Uniprot ID :	O75603

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03388(GCM2 antibody) at dilution of 1:500