

KAT2A antibody

Product Information

Catalog No.:	FNab03389
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	K(lysine) acetyltransferase 2A
Synonyms:	Histone acetyltransferase KAT2A General control of amino acid synthesis protein 5-like 2 Histone acetyltransferase GCN5 (hGCN5) Histone glutaryltransferase KAT2A Histone succinyltransferase KAT2A Lysine acetyltransferase 2A STAF97 KAT2A GCN5 GCN5L2
Observed MW:	94 kDa
Uniprot ID :	Q92830

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab03389(GCN5 antibody) at dilution of 1:800