

NR6A1 antibody

Product Information

Catalog No.:	FNab03390
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	nuclear receptor subfamily 6, group A, member 1
Synonyms:	Nuclear receptor subfamily 6 group A member 1 Germ cell nuclear factor (GCNF, hGCNF) Retinoid receptor-related testis-specific receptor (RTR, hRTR) NR6A1 GCNF
Observed MW:	55 kDa
Uniprot ID :	Q15406

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis tissue slide using FNab03390(NR6A1 Antibody) at dilution of 1:50



IP Result of anti-GCNF (IP:FNab03390, 4ug; Detection:FNab03390 1:800) with mouse testis tissue lysate 4000ug.



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab03390(NR6A1 antibody) at dilution of 1:300