

CSF3R antibody

Product Information

Catalog No.:	FNab03396
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	colony stimulating factor 3 receptor(granulocyte)
Synonyms:	Granulocyte colony-stimulating factor receptor (G-CSF receptor, G-CSF-R) CSF3R GCSFR
Observed MW:	85 kDa
Uniprot ID :	Q99062

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human placenta tissue were subjected to SDS
PAGE followed by western blot with
FNab03396(G-CSFR antibody) at dilution of 1:300