

GDI2 antibody

Product Information

Catalog No.:	FNab03407
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Regulates the GDP/GTP exchange reaction of most Rab proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them.

Immunogen information


Immunogen:	GDP dissociation inhibitor 2
Synonyms:	Rab GDP dissociation inhibitor beta (Rab GDI beta) Guanosine diphosphate dissociation inhibitor 2 (GDI-2) GDI2 RABGDIB
Observed MW:	46-50 kDa
Uniprot ID :	P50395

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab03407(GDI2 antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03407(GDI2 antibody) at dilution of 1:300