

GDI2 antibody

Product Information

Catalog No.:	FNab03408
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0B10
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Regulates the GDP/GTP exchange reaction of most Rab proteins by inhibiting the dissociation of GDP from them, and the subsequent binding of GTP to them.

Immunogen information


Immunogen:	GDP dissociation inhibitor 2
Synonyms:	Rab GDP dissociation inhibitor beta (Rab GDI beta) Guanosine diphosphate dissociation inhibitor 2 (GDI-2) GDI2 RABGDIB
Observed MW:	46 kDa
Uniprot ID :	P50395

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using FNab03408(GDI2 Antibody) at dilution of 1:50



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03408(GDI2 antibody) at dilution of 1:1000