

## ARHGEF25 antibody

### Product Information

Catalog No.:	FNab03413
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	RhoA/RAC/CDC42 exchange factor
Synonyms:	Rho guanine nucleotide exchange factor 25 Guanine nucleotide exchange factor GEFT Rac/Cdc42/Rho exchange factor GEFT RhoA/Rac/Cdc42 guanine nucleotide exchange factor GEFT p63RhoGEF ARHGEF25 GEFT
Observed MW:	64 kDa
Uniprot ID :	Q86VW2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-GEFT (IP:FNab03413, 4ug;  
Detection:FNab03413 1:400) with HEK-293 cells  
lysate 2000ug.



Mouse skeletal muscle tissue were subjected to  
SDS PAGE followed by western blot with  
FNab03413(GEFT antibody) at dilution of 1:1500