

ARHGEF25-Specific antibody

Product Information

Catalog No.:	FNab03416
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	RhoA/RAC/CDC42 exchange factor
Synonyms:	Rho guanine nucleotide exchange factor 25 Guanine nucleotide exchange factor GEFT Rac/Cdc42/Rho exchange factor GEFT RhoA/Rac/Cdc42 guanine nucleotide exchange factor GEFT p63RhoGEF ARHGEF25 GEFT
Observed MW:	64 kDa
Uniprot ID :	Q86VW2

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IP
Recommended dilution:	WB: 1:500-1:3000; IP: 1:500-1:3000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas slide using FNab03416 (GEFT-Specific Antibody) at dilution of 1:100



IP Result of anti-GEFT-Specific (IP:FNab03416, 4ug; Detection:FNab03416 1:1500) with mouse brain tissue lysate 8000ug.



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with FNab03416 (GEFT-Specific antibody) at dilution of 1:1500