

## GSN antibody

### Product Information

Catalog No.:	FNab03417
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	gelsolin(amyloidosis, Finnish type)
Synonyms:	Gelsolin AGEL Actin-depolymerizing factor (ADF) Brevin GSN
Observed MW:	87-93 kDa
Uniprot ID :	P06396

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03417(GSN antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03417(Gelsolin Antibody) at dilution of 1:1000