

GMNN antibody

Product Information

Catalog No.:	FNab03424
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	geminin, DNA replication inhibitor
Synonyms:	Geminin GMNN
Observed MW:	28-29 kDa
Uniprot ID :	O75496

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP, FC, IHC
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using FNab03424(GMNN Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



IP Result of anti-Geminin (IP:FNab03424, 4ug; Detection:FNab03424 1:1000) with HEK-293 cells lysate 4500ug.



human testis tissue were subjected to SDS PAGE followed by western blot with FNab03424(GMNN antibody) at dilution of 1:300