

GOLGB1 antibody

Product Information

Catalog No.:	FNab03453
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	golgin B1, golgi integral membrane protein
Synonyms:	Golgin subfamily B member 1 372 kDa Golgi complex-associated protein (GCP372) Giantin Macrogolgin GOLGB1
Observed MW:	376 kDa
Uniprot ID :	Q14789

Application

Reactivity:	Human
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:20-1:200
Image:	



Immunofluorescent analysis of HepG2 cells using FNab03453(Giantin antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated Goat Anti-Rabbit IgG(H+L)