

GIMAP5 antibody

Product Information

Catalog No.:	FNab03458
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Required for mitochondrial integrity and T-cell survival. May contribute to T-cell quiescence(By similarity).

Immunogen information


Immunogen:	GTPase, IMAP family member 5
Synonyms:	GTPase IMAP family member 5 Immune-associated nucleotide-binding protein 5 Immunity-associated nucleotide 4-like 1 protein Immunity-associated nucleotide 5 protein (IAN-5, hIAN5) Immunity-associated protein 3 GIMAP5 IAN4L1 IAN5 IMAP3
Observed MW:	45-50 kDa
Uniprot ID :	Q96F15

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lymphoma using FNab03458(GIMAP5 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab03458(GIMAP5 antibody) at dilution of 1:500