

GJB6 antibody

Product Information

Catalog No.:	FNab03468
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Immunogen information

Immunogen:	gap junction protein, beta 6, 30kDa
Synonyms:	Gap junction beta-6 protein Connexin-30 (Cx30) GJB6
Observed MW:	None
Uniprot ID :	O95452

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human oesophagus tissue slide using FNab03468(GJB6 Antibody) at dilution of 1:50