

GK3 antibody

Product Information

Catalog No.:	FNab03470
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Key enzyme in the regulation of glycerol uptake and metabolism.

Immunogen information


Immunogen:	glycerol kinase 3 pseudogene
Synonyms:	Glycerol kinase 3 (GK 3, Glycerokinase 3) ATP:glycerol 3-phosphotransferase 3 Glycerol kinase 3 pseudogene Glycerol kinase, testis specific 1 GK3 GK3P GKP3 GKTB
Observed MW:	61 kDa
Uniprot ID :	Q14409

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human skin cancer using FNab03470(GK3P antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03470(GK3P antibody) at dilution of 1:500