

## GKN1 antibody

### Product Information

Catalog No.:	FNab03473
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6F7
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Has mitogenic activity and may be involved in maintaining the integrity of the gastric mucosal epithelium.

### Immunogen information


Immunogen:	gastrokin 1
Synonyms:	Gastrokin-1 18 kDa antrum mucosa protein (AMP-18) Protein CA11 GKN1 AMP18 CA11
Observed MW:	22 kDa
Uniprot ID :	Q9NS71

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:5000-1:50000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human stomach tissue slide using FNab03473( GKN1 Antibody) at dilution of 1:200



human stomach tissue were subjected to SDS PAGE followed by western blot with FNab03473( GKN1 Antibody) at dilution of 1:20000