

GLI3-Specific antibody

Product Information


Catalog No.:	FNab03478
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	GLI family zinc finger 3
Synonyms:	Transcriptional activator GLI3 GLI3 form of 190 kDa (GLI3-190) GLI3 full-length protein (GLI3FL) Transcriptional repressor GLI3R Alternative names: GLI3 C-terminally truncated form GLI3 form of 83 kDa (GLI3-83) GLI3
Observed MW:	230 kDa
Uniprot ID :	P10071

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:200-1:1000
Image:	



human placenta tissue were subjected to SDS
PAGE followed by western blot with
FNab03478(GLI3-Specific antibody) at dilution of
1:400