

QARS1 antibody

Product Information

Catalog No.:	FNab03482
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	glutaminyt-tRNA synthetase
Synonyms:	Glutamine--tRNA ligase Glutaminyt-tRNA synthetase (GlnRS) QARS1 QARS
Observed MW:	88 kDa
Uniprot ID :	P47897

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:3000; IHC: 1:20-1:200; IF: 1:10-1:100
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab03482(QARS antibody) at dilution of 1:50



PC-3 cells were subjected to SDS PAGE followed by western blot with FNab03482(QARS antibody) at dilution of 1:500