

COLGALT1 antibody

Product Information

Catalog No.:	FNab03493
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	glycosyltransferase 25 domain containing 1
Synonyms:	Procollagen galactosyltransferase 1 Collagen beta(1-O)galactosyltransferase 1 (ColGalT 1) Glycosyltransferase 25 family member 1 Hydroxylysine galactosyltransferase 1 COLGALT1 GLT25D1
Observed MW:	72 kDa
Uniprot ID :	Q8NBJ5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03493(GLT25D1 Antibody) at dilution of 1:1000