

COLGALT2 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab03494 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information

| | |
|--------------|--|
| Immunogen: | glycosyltransferase 25 domain containing 2 |
| Synonyms: | Procollagen galactosyltransferase 2 Collagen beta(1-O)galactosyltransferase 2 (ColGalT 2) Glycosyltransferase 25 family member 2 Hydroxylysine galactosyltransferase 2 COLGALT2 C1orf17 GLT25D2 KIAA0584 |
| Observed MW: | 65 kDa |
| Uniprot ID : | Q8IYK4 |

Application

| | |
|--------------------------|--------------|
| Reactivity: | Human, Mouse |
| Tested Application: | ELISA, WB |
| Recommended dilution:WB: | 1:500-1:2000 |
| Image: | |



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab03494(GLT25D2 Antibody) at dilution of 1:600