

GCG antibody

Product Information

Catalog No.:	FNab03496
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	glucagon
Synonyms:	Pro-glucagon Glicentin Glicentin-related polypeptide (GRPP) Oxyntomodulin (OXM, OXY) Glucagon Glucagon-like peptide 1 (GLP-1) Alternative names: Incretin hormone Glucagon-like peptide 1(7-37) (GLP-1(7-37)) Glucagon-like peptide 1(7-36) (GLP-1(7-36)) Glucagon-like peptide 2 (GLP-2) GCG
Observed MW:	21 kDa
Uniprot ID :	P01275

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemical of paraffin-embedded human pancreas using FNab03496(GCG antibody) at dilution of 1:200



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab03496(Glucagon Antibody) at dilution of 1:500