

GLUL antibody

Product Information

Catalog No.:	FNab03505
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	glutamate-ammonia ligase(glutamine synthetase)
Synonyms:	Glutamine synthetase (GS) Glutamate--ammonia ligase Palmitoyltransferase GLUL GLUL GLNS
Observed MW:	44 kDa
Uniprot ID :	P15104

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, FC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:100-1:400
Image:	



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using FNab03505 (Glutamine synthetase Antibody) at dilution of 1:200 heat mediated antigen retrieved with Tris-EDTA buffer (pH9).



rat liver tissue were subjected to SDS PAGE followed by western blot with FNab03505 (Glutamine synthetase Antibody) at dilution of 1:600