

GPC4 antibody

Product Information

Catalog No.:	FNab03514
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	glypican 4
Synonyms:	Glypican-4 K-glypican Secreted glypican-4 GPC4
Observed MW:	62 kDa
Uniprot ID :	O75487

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab03514(Glypican 4 antibody) at dilution of 1:200