

GMFG antibody

Product Information

Catalog No.:	FNab03524
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	glia maturation factor, gamma
Synonyms:	Glia maturation factor gamma (GMF-gamma) GMFG
Observed MW:	17 kDa
Uniprot ID :	O60234

Application


Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human lung cancer tissue slide using FNab03524(GMFG Antibody) at dilution of 1:50



IP Result of anti-GMFG (IP:FNab03524, 4ug; Detection:FNab03524 1:500) with mouse spleen tissue lysate 4000ug.



rat brain tissue were subjected to SDS PAGE followed by western blot with FNab03524(GMFG Antibody) at dilution of 1:600