

## GOSR2 antibody

### Product Information

Catalog No.:	FNab03565
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	golgi SNAP receptor complex member 2
Synonyms:	Golgi SNAP receptor complex member 2 27 kDa Golgi SNARE protein Membrin GOSR2 GS27
Observed MW:	25 kDa
Uniprot ID :	O14653

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



JurKat cells were subjected to SDS PAGE followed by western blot with FNab03565(GOSR2/Membrin antibody) at dilution of 1:1000