

GOSR2 antibody

Product Information

Catalog No.:	FNab03566
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8G6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


This gene encodes a trafficking membrane protein which transports proteins among the medial- and trans-Golgi compartments. Due to its chromosomal location and trafficking function, this gene may be involved in familial essential hypertension.

Immunogen information


Immunogen:	golgi SNAP receptor complex member 2
Synonyms:	Golgi SNAP receptor complex member 2 27 kDa Golgi SNARE protein Membrin GOSR2 GS27
Observed MW:	25 kDa
Uniprot ID :	O14653

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver slide using FNab03566(GOSR2 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03566(GOSR2 antibody) at dilution of 1:1000