

## GPAA1 antibody

### Product Information

Catalog No.:	FNab03574
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Essential for GPI-anchoring of precursor proteins but not for GPI synthesis. Acts before or during formation of the carbonyl intermediate.

### Immunogen information


Immunogen:	glycosylphosphatidylinositol anchor attachment protein 1 homolog(yeast)
Synonyms:	Glycosylphosphatidylinositol anchor attachment 1 protein (GPI anchor attachment protein 1) GAA1 protein homolog (hGAA1) GPAA1 GAA1
Observed MW:	60-68 kDa
Uniprot ID :	O43292

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab03574(GPAA1 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab03574(GPAA1 antibody) at dilution of 1:400