

P2RY14 antibody

Product Information

Catalog No.:	FNab03588
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	purinergic receptor P2Y, G-protein coupled, 14
Synonyms:	P2Y purinoceptor 14 (P2Y14) G-protein coupled receptor 105 UDP-glucose receptor P2RY14 GPR105 KIAA0001
Observed MW:	None
Uniprot ID :	Q15391

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab03588(GPR105 antibody) at dilution of 1:50