

GPR142 antibody

Product Information

Catalog No.:	FNab03597
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Orphan receptor.

Immunogen information

Immunogen:	G protein-coupled receptor 142
Synonyms:	Probable G-protein coupled receptor 142 G-protein coupled receptor PGR2 GPR142 PGR2
Observed MW:	51-55 kDa
Uniprot ID :	Q7Z601

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab03597(GPR142 antibody) at dilution of 1:500