

## GPR183 antibody

### Product Information

Catalog No.:	FNab03602
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	G protein-coupled receptor 183
Synonyms:	G-protein coupled receptor 183 Epstein-Barr virus-induced G-protein coupled receptor 2 (EBI2, EBV-induced G-protein coupled receptor 2, hEBI2) GPR183 EBI2
Observed MW:	55-60 kDa
Uniprot ID :	P32249

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-GPR183/EBI2 (IP:FNab03602, 4ug; Detection:FNab03602 1:500) with mouse lung tissue lysate 4000ug.



JurKat cells were subjected to SDS PAGE followed by western blot with FNab03602(GPR183/EBI2 antibody) at dilution of 1:1000