

GPR3 antibody

Product Information

Catalog No.:	FNab03604
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Orphan receptor with constitutive G(s) signaling activity that activate cyclic AMP. Has a potential role in modulating a number of brain functions, including behavioral responses to stress(By similarity), amyloid-beta peptide generation in neurons and neurite outgrowth(By similarity). Maintains also meiotic arrest in oocytes(By similarity).

Immunogen information

Immunogen:	G protein-coupled receptor 3
Synonyms:	G-protein coupled receptor 3 ACCA orphan receptor GPR3 ACCA
Observed MW:	None
Uniprot ID :	P46089

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human kidney tissue slide using FNab03604(GPR3 antibody) at dilution of 1:50