

GPR37 antibody

Product Information

Catalog No.:	FNab03606
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	G protein-coupled receptor 37(endothelin receptor type B-like)
Synonyms:	Prosaposin receptor GPR37 Endothelin B receptor-like protein 1 (ETBR-LP-1) G-protein coupled receptor 37 Parkin-associated endothelin receptor-like receptor (PAELR) GPR37
Observed MW:	50-55 kDa
Uniprot ID :	O15354

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03606(GPR37/Pael-R antibody) at dilution of 1:300