

GPR50 antibody

Product Information

Catalog No.:	FNab03608
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Does not bind melatonin.

Immunogen information


Immunogen:	G protein-coupled receptor 50
Synonyms:	Melatonin-related receptor G protein-coupled receptor 50 H9 C-terminal domain GPR50
Observed MW:	67 kDa
Uniprot ID :	Q13585

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IF,IP, FC
Recommended dilution:	IP: 1:200-1:2000; IF: 1:20-1:200
Image:	



Immunofluorescent analysis of SH-SY5Y cells using FNab03608(GPR50Antibody) at dilution of 1:100



IP Result of anti-GPR50 (IP:FNab03608, 4ug; Detection:FNab03608 1:500) with Hela cell lysate 2000ug.