

TAAR2 antibody

Product Information

Catalog No.:	FNab03609
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	trace amine associated receptor 2
Synonyms:	Trace amine-associated receptor 2 (TaR-2, Trace amine receptor 2) G-protein coupled receptor 58 TAAR2 GPR58
Observed MW:	42 kDa
Uniprot ID :	Q9P1P5

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Rat brain tissue were subjected to SDS PAGE followed by western blot with FNab03609(GPR58/TAAR2 antibody) at dilution of 1:1000