

GPR83 antibody

Product Information

Catalog No.:	FNab03611
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Orphan receptor. Could be a neuropeptide Y receptor.

Immunogen information

Immunogen:	G protein-coupled receptor 83
Synonyms:	G-protein coupled receptor 83 G-protein coupled receptor 72 GPR83 GPR72 JP05 KIAA1540
Observed MW:	48-55 kDa
Uniprot ID :	Q9NYM4

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03611(GPR83 antibody) at dilution of 1:500