

GPRIN1 antibody

Product Information

Catalog No.:	FNab03615
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in neurite outgrowth.

Immunogen information


Immunogen:	G protein regulated inducer of neurite outgrowth 1
Synonyms:	G protein-regulated inducer of neurite outgrowth 1 (GRIN1) GPRIN1 KIAA1893
Observed MW:	100-110 kDa
Uniprot ID :	Q7Z2K8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, IP, WB
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:1000
Image:	



IP Result of anti-GPRIN1 (IP:FNab03615, 3ug;
Detection:FNab03615 1:500) with mouse brain
tissue lysate 3600ug.



mouse brain tissue were subjected to SDS PAGE
followed by western blot with
FNab03615(GPRIN1 antibody) at dilution of
1:1000