

## GR repeat antibody

### Product Information

Catalog No.:	FNab03626
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	None
Synonyms:	None
Observed MW:	18 kDa/36 kDa
Uniprot ID :	None

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Recombinant protein were subjected to SDS  
PAGE followed by western blot with  
FNab03626(GR repeat antibody) at dilution of  
1:1000