

## GRN antibody

### Product Information

|               |                                                                                                                |
|---------------|----------------------------------------------------------------------------------------------------------------|
| Catalog No.:  | FNab03631                                                                                                      |
| Size:         | 100µg                                                                                                          |
| Form:         | liquid                                                                                                         |
| Purification: | Immunogen affinity purified                                                                                    |
| Purity:       | ≥95% as determined by SDS-PAGE                                                                                 |
| Host:         | Rabbit                                                                                                         |
| Clonality:    | polyclonal                                                                                                     |
| Clone ID:     | None                                                                                                           |
| IsoType:      | IgG                                                                                                            |
| Storage:      | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

### Background

None


### Immunogen information

|              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen:   | granulin                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Synonyms:    | Progranulin (PGRN) Acrogranin Epithelin precursor Glycoprotein of 88 Kda (GP88, Glycoprotein 88) Granulin precursor PC cell-derived growth factor (PCDGF) Proepithelin (PEPI) Paragranulin Granulin-1<br>Alternative names: Granulin G Granulin-2 Alternative names: Granulin F Granulin-3 Alternative names: Epithelin-2 Granulin B Granulin-4<br>Alternative names: Epithelin-1 Granulin A Granulin-5 Alternative names: Granulin C Granulin-6 Alternative names: Granulin D Granulin-7 Alternative names: Granulin E GRN |
| Observed MW: | 78 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Uniprot ID : | P28799                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |


### Application

|                       |                                   |
|-----------------------|-----------------------------------|
| Reactivity:           | Human                             |
| Tested Application:   | ELISA, WB, IHC                    |
| Recommended dilution: | WB: 1:500-1:1000; IHC: 1:20-1:200 |

Image:



Immunohistochemistry of paraffin-embedded human prostate cancer slide using FNab03631 (PCDGF, GRN Antibody) at dilution of 1:50



Hela cells were subjected to SDS PAGE followed by western blot with FNab03631 (PCDGF, GRN antibody) at dilution of 1:1000