

GRIPAP1 antibody

Product Information

Catalog No.:	FNab03651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	GRIP1 associated protein 1
Synonyms:	GRIP1-associated protein 1 (GRASP-1) GRASP-1 C-terminal chain Alternative names: 30kDa C-terminus form GRIPAP1 KIAA1167
Observed MW:	97 kDa
Uniprot ID :	Q4V328

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab03651 (GRIPAP1 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03651 (GRIPAP1 antibody) at dilution of 1:500