

GOSR1 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab03668 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


None

Immunogen information


| | |
|--------------|--|
| Immunogen: | golgi SNAP receptor complex member 1 |
| Synonyms: | Golgi SNAP receptor complex member 1 28 kDa Golgi SNARE protein 28 kDa cis-Golgi SNARE p28 (GOS-28) GOSR1 GS28 |
| Observed MW: | 28 kDa |
| Uniprot ID : | O95249 |

Application

| | |
|-----------------------|-----------------------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB, IHC, FC |
| Recommended dilution: | WB: 1:500-1:2000; IHC: 1:20-1:200 |
| Image: | |



Immunohistochemistry of paraffin-embedded human stomach cancer using FNab03668(GS28 antibody) at dilution of 1:50



human liver tissue were subjected to SDS PAGE followed by western blot with FNab03668(GS28 antibody) at dilution of 1:500