

GSTT1 antibody

Product Information

Catalog No.:	FNab03700
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Acts on 1,2-epoxy-3-(4-nitrophenoxy)propane, phenethylisothiocyanate 4-nitrobenzyl chloride and 4-nitrophenethyl bromide. Displays glutathione peroxidase activity with cumene hydroperoxide.

Immunogen information


Immunogen:	glutathione S-transferase theta 1
Synonyms:	Glutathione S-transferase theta-1 GST class-theta-1 Glutathione transferase T1-1 GSTT1
Observed MW:	28 kDa, 30 kDa
Uniprot ID :	P30711

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab03700(GSTT1 antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03700(GSTT1 antibody) at dilution of 1:500