

## OLA1 antibody

### Product Information

Catalog No.:	FNab03724
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	Obg-like ATPase 1
Synonyms:	Obg-like ATPase 1 DNA damage-regulated overexpressed in cancer 45 (DOC45) GTP-binding protein 9 OLA1 GTPBP9
Observed MW:	45 kDa
Uniprot ID :	Q9NTK5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab03724(OLA1 antibody) at dilution of 1:50



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab03724(GTPBP9 antibody) at dilution of 1:600