

GUCY1A1 antibody

Product Information

Catalog No.:	FNab03726
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	guanylate cyclase 1, soluble, alpha 3
Synonyms:	Guanylate cyclase soluble subunit alpha-1 (GCS-alpha-1) Guanylate cyclase soluble subunit alpha-3 (GCS-alpha-3) Soluble guanylate cyclase large subunit GUCY1A1 GUC1A3 GUCSA3 GUCY1A3
Observed MW:	77 kDa
Uniprot ID :	Q02108

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab03726(GUCY1A3 Antibody) at dilution of 1:50



human brain tissue were subjected to SDS PAGE followed by western blot with FNab03726(GUCY1A3 antibody) at dilution of 1:400