

GYG2 antibody

Product Information

Catalog No.:	FNab03735
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Self-glucosylates, via an inter-subunit mechanism, to form an oligosaccharide primer that serves as substrate for glycogen synthase.

Immunogen information

Immunogen:	glycogenin 2
Synonyms:	Glycogenin-2 (GN-2, GN2) GYG2
Observed MW:	45 kDa, 66-70 kDa
Uniprot ID :	O15488

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



L02 cells were subjected to SDS PAGE followed by western blot with FNab03735(GYG2 antibody) at dilution of 1:500