

HACL1 antibody

Product Information

Catalog No.:	FNab03746
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes a carbon-carbon cleavage reaction; cleaves a 2-hydroxy-3-methylacyl-CoA into formyl-CoA and a 2-methyl-branched fatty aldehyde.

Immunogen information


Immunogen:	2-hydroxyacyl-CoA lyase 1
Synonyms:	2-hydroxyacyl-CoA lyase 1 2-hydroxyphytanoyl-CoA lyase (2-HPCL) Phytanoyl-CoA 2-hydroxylase 2 HACL1 HPCL HPCL2 PHYH2
Observed MW:	64 kDa
Uniprot ID :	Q9UJ83

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab03746(HAACL1 antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab03746(HAACL1 antibody) at dilution of 1:500