

HAPLN4 antibody

Product Information

Catalog No.:	FNab03754
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to hyaluronic acid and may be involved in formation of the extracellular matrix.

Immunogen information


Immunogen:	hyaluronan and proteoglycan link protein 4
Synonyms:	Hyaluronan and proteoglycan link protein 4 Brain link protein 2 HAPLN4 BRAL2 KIAA1926
Observed MW:	43 kDa
Uniprot ID :	Q86UW8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab03754(HAPLN4 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab03754(HAPLN4 antibody) at dilution of 1:1200