

HARS1 antibody

Product Information

Catalog No.:	FNab03758
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Cytoplasmic histidine--tRNA ligase(Probable). Plays a role in axon guidance.

Immunogen information

Immunogen:	histidyl-tRNA synthetase
Synonyms:	Histidine--tRNA ligase, cytoplasmic Histidyl-tRNA synthetase (HisRS) HARS1 HARS HRS
Observed MW:	50 kDa
Uniprot ID :	P12081

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab03758(HARS antibody) at dilution of 1:500