

## MNX1 antibody

### Product Information

Catalog No.:	FNab03766
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	motor neuron and pancreas homeobox 1
Synonyms:	Motor neuron and pancreas homeobox protein 1 Homeobox protein HB9 MNX1 HLXB9
Observed MW:	48 kDa
Uniprot ID :	P50219

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:1000
Image:	



mouse spleen tissue were subjected to SDS PAGE followed by western blot with FNab03766( HB9 Antibody) at dilution of 1:600