

HBG1-Specific antibody

Product Information

Catalog No.:	FNab03772
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0H6
IsoType:	IgG2a
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	hemoglobin, gamma A
Synonyms:	Hemoglobin subunit gamma-1 Gamma-1-globin Hb F Agamma Hemoglobin gamma-1 chain Hemoglobin gamma-A chain HBG1
Observed MW:	16 kDa
Uniprot ID :	P69891

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab03772 (HBG1-Specific antibody) at dilution of 1:200