

HBQ1 antibody

Product Information

Catalog No.:	FNab03775
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	hemoglobin, theta 1
Synonyms:	Hemoglobin subunit theta-1 Hemoglobin theta-1 chain Theta-1-globin HBQ1
Observed MW:	None
Uniprot ID :	P09105

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human hepatocirrhosis using FNab03775(HBQ1 antibody) at dilution of 1:50