

## HCFC2 antibody

### Product Information

Catalog No.:	FNab03783
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	host cell factor C2
Synonyms:	Host cell factor 2 (HCF-2) C2 factor HCFC2
Observed MW:	87 kDa
Uniprot ID :	Q9Y5Z7

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab03783(HCFC2 antibody) at dilution of 1:300